Anclasil oclu

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E

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Editing date: 13.09.2010

Instructions for use

**Bite registration material, silicone based, addition curing
fast setting – extreme high final hardness**

Fig. 1

1. Mixing and Dosing

Extrude material from the cartridge with Automix2-system gun (fig. 1). Place cartridge into mixing gun and remove closure cap by rotation. Eject a small amount of material (bleeding) until equal flow is being extruded from both orifices (fig. 2). Insert mixing cannula (green) into guide grooves on cartridge and cannula and rotate anti-clockwise to a stop (fig. 3). Store cartridge with mixing cannula attached for subsequent use.


Fig. 2

2. Instruction for use

Extrude material from cartridge with mixing cannula or with attached contouring tip onto teeth (fig. 4). At terminal occlusion application from vestibular is possible. Let registration harden in terminal occlusion and remove from the mouth after 1 min.

Fig. 3

3. Desinfection


After removal from the mouth rinse registrate under lukewarm water. Disinfection for 15 min. in 2 % glutaraldehyde.

Fig. 4

Important working hints

- Remove any retraction astringents or mouth wash solutions prior to the registration – could affect setting reaction.
- Do not intermix with condensation curing silicones.
- Avoid staining of clothing. Hardened impression materials cannot be removed.
- Remove all residues of registration material from the patient's mouth.
- For reducing the pressure needed during injection, we recommend using the green mixing cannula.

Indications for use:

- Bite registration with crown and bridge impressions and for orthodontic indications, key for intraoral registration

Technical Data:

- DIN 13903
- **Mixed volume:** 50 ml (cartridge)
- **Mixing ratio:** 1:1
- **Colour code:** base: green catalyst: yellow
- **Mixing time:** omitted (Automix2-system)
- **Working time:** approx. 30 sec.*
- **Setting time:** approx. 1 min. 30 sec.*
- **Time in mouth:** approx. 1 min.
- **Hardness Shore D:** 45
- **Linear dimensional change:** < 0.2 %
- **Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.

Ordering information:

Anclasil ocu A-040014
2 x 50 ml cartridges
12 mixing cannulas
12 contouring tips

Anclasil ocu A-040015
mixing cannulas
Bag with 25 mixing cannulas

Anclasil ocu A-040016
contouring tips
Bag with 25 contouring tips

Ancladén

Ancladén, S.L. | Rubinstein, 2 Entlo 1ª | 08022 Barcelona
93 417 29 84 | ancladen@ancladen.com | www.ancladen.com