Anclacem temp

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Anclacem temp

Instructions for use

Temporary luting cement, dual-hardening, for crowns and bridges, eugenol-free – composite with triclosan


Fig. 1

1. Cleaning the preparation

Before cementing, clean the preparations with water. If the preparation has been disinfected beforehand, remove all traces of disinfectant from the surface. Dry the preparation area and dry the temporary restoration (Fig. 1).


Fig. 2

2. Mixing and applying

Dispensing is performed with the mini-mix system. Insert plunger into the body of the cartridge. Remove the cap from the cartridge by turning. Before fitting the mixing cannula, dispense a small quantity of material until material emerges from the two openings uniformly (Fig. 2). When attaching the mixing cannula observe the guides on the cannula and cartridge. Lock in place by turning in the opposite direction (Fig. 3). Before each application dispense and discard a small quantity. Dispensing can now be performed to suit requirements. Dispense material with uniform pressure. After use, leave the mixing cannula on the cartridge until the next application.

Fig. 3

3. Inserting the temporary restoration

Apply a thin layer of **Anclacem temp** to the interior of the temporary restorations (Fig. 4). Directly after that, position the temporary restorations on the preparations by applying slight pressure (Fig. 5), then press them on firmly. Remove any excess material with a suitable instrument (Fig. 6).


Fig. 4

4. Working time / light curing

The working time at room temperature is approx. 90 sec. Curing at body temperature takes approx. 5 min. and can be delayed outside the mouth. **Setting can be assisted and speeded up at any time by light curing.** For this purpose we recommend an exposure time of 20 sec. per surface (Fig. 7). If the temporary restorations are very opaque, extend light curing to 30 - 40 sec.


Fig. 5

Fig. 6

Important working hints

- **Anclacem temp** is compatible with permanent cements.
- **Anclacem temp** does not have a delaying effect on the curing of composite cement or composite restorative material.
- In the case of composite-based and build-up restorations isolate these areas by applying a thin layer of glycerine or Vaseline in order to stop all adhesion of **Anclacem temp**.
- Before using **Anclacem temp** do not apply any methacrylate-based bonder, dentine sealant or desensitizer.
- **Anclacem temp** does not contain eugenol.
- In very rare cases **Anclacem temp** may cause sensitization. If such reactions should occur, discontinue using **Anclacem temp**.
- Manual removal and mixing of the components is not recommended.
- Do not allow individual components to come into contact with skin or mucous membrane.
- Do not leave any residues of material in the mouth.


Fig. 7

Dual-hardening system

light curing:

self curing:

Applications:

- Temporary luting of crowns and bridges
- Temporary luting of permanent restorations

Technical data:

- **Mixed volume:** 5 ml (mini-mix)
- **Mixing ratio:** 4:1
- **Colour code:** Base: transparent
Catalyst: transparent
- **Mixing time:** omitted (mini-mix)
- **Working time:** approx. 90 sec.
- **Setting time in the mouth:** light curing: 20 sec. per surface
self curing: 5 min.
- **Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.

Ordering information:

Anclacem temp A-050005
cartridge of 5 ml mini-mix 4:1
8 mixing cannulas 4:1,
brown

Mixing cannulas 4:1
25 pcs. **A-050002**

Ancladén

Ancladén, S.L. | Rubinstein, 2 Entlo 1ª | 08022 Barcelona
93 417 29 84 | ancladen@ancladen.com | www.ancladen.com